

Papiers-filtres standard disques à plat

vitesse	moyenne	rapide	très rapide
poids	73 g/m ²	64 g/m ²	73 g/m ²
résistance	1,95 kg/cm ²	1,45 kg/cm ²	1,25 kg/cm ²
épaisseur papier	160 µm	150 µm	170 µm
temps de filtration	88 s	50 s	22 s
porosité	5 à 13 µm	10 à 20 µm	17 à 30 µm
Ø 20 mm		FL1231	
Ø 25 mm	FL1402	FL1232	FL1411
Ø 32 mm		FL1233	
Ø 35 mm	FL1403		
Ø 37 mm	FL1404		
Ø 40 mm	FL1405		
Ø 42,5 mm	FL1380	FL1234	-
Ø 45 mm	FL1381	FL1210	FL1280
Ø 47 mm	FL1382	FL1211	FL1412
Ø 50 mm	FL1383	-	-
Ø 55 mm	FL1384	FL1212	FL1281
Ø 60 mm			FL1413
Ø 70 mm	FL1385	FL1213	FL1282
Ø 80 mm		FL1235	FL1414
Ø 85 mm	FL1386	FL1236	-
Ø 90 mm	FL1387	FL1214	FL1283
Ø 100 mm	FL1388	FL1215	FL1284
Ø 110 mm	FL1389	FL1216	FL1285
Ø 120 mm	FL1406	FL1237	FL1415
Ø 125 mm	FL1390	FL1217	FL1286
Ø 130 mm	FL1391	FL1218	FL1287
Ø 150 mm	FL1392	FL1219	FL1288
Ø 180 mm		FL1238	FL1416
Ø 185 mm	FL1393	FL1220	FL1289
Ø 190 mm	FL1394	FL1221	FL1290
Ø 200 mm	-	FL1222	FL1291
Ø 210 mm	FL1395	FL1223	FL1292
Ø 230 mm	FL1407	FL1239	
Ø 240 mm	FL1396	FL1224	FL1293
Ø 250 mm	FL1408	FL1225	FL1294
Ø 270 mm	FL1397	FL1254	FL1295
Ø 300 mm		FL1255	
Ø 305 mm			FL1417
Ø 320 mm	FL1398	FL1226	FL1418
Ø 330 mm	-	FL1227	FL1296
Ø 360 mm	FL1409		
Ø 400 mm	FL1399	FL1228	FL1297
Ø 500 mm	FL1400	FL1230	FL1298
Ø 600 mm	FL1410		
Ø 620 mm			FL1419
Ø 1200 mm	FL1401	-	-



- papiers-filtres lisses en cellulose pure
- poids selon DIN 53104
- temps de filtration selon DIN 53137
- la boîte de 100 filtres